

Luggie SUPER

The Luggie scooter is your ideal carry-on scooter. Convenient for short distances, foldable under a minute and super easy to carry-on like a travelling suitcase on a boat, train or airplane (including its battery!). Handy to negotiate sharp turns in small rooms, restaurants or shops. Complete with adjustable and removable armrests.



Completely foldable and easy to transport in your car

1. Fold,
2. extend steering column
3. go...



users weight
164 KG

- Length : 989 mm
- Width : 525 mm
- Height : 860 mm
- Distance ground - footplate : 120 mm
- Battery : 10,5 Ah
- Driving range : 18 km**
- Weight: 30,5 kg
- Speed : 6 km/h*
- Users weight : 164 kg
- Color : blue



Model shown only for illustration purposes. Certain features may be optional.

www.vermeiren.com - info@vermeiren.com