






* Use:



WEICONLOCK AN 302-43

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

Sheet 1/6

PLEGO				Plan No.
				E00052670/E.1
Scale: 1:10	Drafted	Verified	Confirmed	
Date	27/09/2023	17/02/2025	17/02/2025	
Initials	BK	WM	MK	





NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
 PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
 WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
 EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
 FROM VERMEIREN GROUP.


Sheet 3/6

PLEGO			
Scale: 1:10	Drafted	Verified	Confirmed
Date	27/09/2023	17/02/2025	17/02/2025
Initials	BK	AS	WM


Plan No.

E00052670/E.1






B (1:5)



C (1:5)





* Use:
WEICONLOCK AN 302-43

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

Sheet 4/6

PLEGO				---
				Plan No.
				E00052670/E.1
Scale: 1:12	Drafted	Verified	Confirmed	
Date	27/09/2023	17/02/2025	17/02/2025	
Initials	BK	AS	WM	



* Use:
WEICONLOCK AN 302-43

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

Sheet 5/6

PLEGO				Plan No.
				E00052670/E.1
Scale: 1:12	Drafted	Verified	Confirmed	
Date	27/09/2023	17/02/2025	17/02/2025	
Initials	BK	AS	WM	

No.	Plan	Article	Qty	24	WASHER M6	2994414	16	40	A00056458	2997886	1		
1	E00063648	---	1	25	WASHER M10 DIN125	2998106	4	41	A00056457	2997884	1		
2	E00062349	---	1*					42	E00052668	---	1		
3	E00062350	---	1	26	M8x25	1990501	2	43	STAR KNOB M6x35 VGF	2990215	1		
4	E00062352	---	1	27	NUT M8	2995398	2						
5	E00056068	---	1	28	M6X25 DIN912	2996887	8	44	X001473	1002994	10		
6	E00063644	---	1	29	X000862	1006463	8						
7	E00062355	---	1	30	NUT DIN985 M6	2994063	4						
8	P00053716	26771	2	31	M6x45	2994086	4	46	A00056460	2997888**	1		
9	P00062160	2677665	1	32	P00053220	26742B8 26742A9 26742B2 26742A7 26742D5 2674221	8	47	A00056459	2997885	1		
10	P00061760	2677765	1					48	M5x25	2994083	2		
11	P00053715	2678965	1					49	NUT M5	2994006	2		
12	1920605	1920605/1920606 **	1					50	M6X10	1904683	4		
								51	E00054750	---**	1		
13	XC01046	1919584	2	33	M4X10 DIN7991	2998006	8	52	E00065217/ E00066678	---**/ ---**	1		
14	M6x16	2994067	4	34	X000106	G003218	4	53	P00066352	T001061	2		
15	NUT M6	2994063	4	35	M8 DIN125	2994415	8						
16	1919238	1919238	1	36	M8X60 ISO7380	2995491	2	54	P00066353	T001060	2		
17	1919240	1919240	1	37	M5x12	2994074	4	55	A00066354	1920500**	1		
18	2993352 ELEMENT 1 /2993352 ELEMENT 2	2993352	1	38	P00054831	26779	2	56	E00066356	---**	1		
				39	A00056456	2997887	1	57	E00066564	---**	1		
								58	E00066475	---**	1		
19	P00061690	T001215 T001214 T001212 T001216 T001209 T001208	2					59	A00066482	2998031**	1		
				60	P00066485	1920046**	2						
				61	NUT M8	1996213	2						
				62**	B31 CRUTCH- HOLDER PLAST.25MM	1906616	1						
20	P00053212	2997848	2					Sheet 6/6					---
21	M8 DIN125B	1001556	2					PLEGO					Plan No.
22	M6x20	2994004	4										E00052670/E.1
23	P00053209	2997847	4										

* 1 pair contains 2 pcs -
check AS system

** Optional equipment

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

Scale: 1:10

Date

Initials

Drafted

27/09/2023

BK

Verified

17/02/2025

AS

Confirmed

17/02/2025

WM