


EC_COM1 EC_COM3



FOR PLEGO
EC_MOT1 - M2
EC_MOT2 - M1

#	Reference	Art number	Description	Qty
1	EC_BAT1	1918833	EC BAT BATTERY HOLDER 24V/12AH WITH KEY VERSO	1
2	EC_COM1	1919238	EC COM PG POWER MODULE NVR2	1
3	EC_COM2	1919239	EC COM PG JOYSTICK NVR2	1
4	EC_MOT1	3994669	EC MOT MOTOR 250W L LINUX BL	1
5	EC_MOT2	3994670	EC MOT MOTOR 250W R LINUX BL	1
6	EC_ZIP1	2993110	EC ZIP BINDING BAND NYLON L=300MM	1
7	EC_CON3,EC_CON4	2998595	EC CON CE5	2
8	EC_A1	2998576	EC A PLEGO VERSO MAIN PWR CABLE	1
9	EC_COM3	1919240	EC COM PG SEALING KIT NVR2	1

BK BLACK
RD RED
LBL LIGHT BLUE
BL BLUE
LGR LIGHT GREEN
GR GREEN
VI VILLET
OR ORANGE
YE YELLOW
WH WHITE
LRD LIGHT RED
PK PINK
GY GRAY
/ AND
I OR

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP

BOM EC PLEGO VERSO DRIVE TRAIN

Id: 1/1

Size: A4	Drafted	Verified	Comformed
Date:	2025-02-21	2025-02-21	2025-02-21
Initials:	FM	FM	TK

2998597/B

Plan No.

2998597/B

