




No.	Plan	Article	Qty
1	A00052508	2998018	1
2	1919239	1919239	1
3	M5x12 Black	2994074	2
4	P00052651	2675865	1
5	P00066148	2998008	1
6	P1996154	1996154	1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

2997983

SECONDARY CONTROLLER PLEGO

Plan No.

E00054750/B1.1

Scale:	1:2	Drafted	Verified	Confirmed
Date	08.05.2023	17/02/2025	17/02/2025	
Initials	BK	WM	MK	

