



#	Reference	Art number	Description	Qty
1	EC_CON1,EC_CON2	2995839	EC CON CABLE SHOE 6.3X0,8 FEMALE	2
2	EC_HST1,EC_HST2	T001105	EC HST Heat Shrink Tubing 12	0,03M
3	EC_HST3	2995867	EC HST Mesh for cable 8	0,45M
4	EC_WIRE1	2995864	EC WIRE 1,3mm2 BK	0,5M
5	EC_WIRE2	2995865	EC WIRE 1,3mm2 RD	0,5M

BK	BLACK
RD	RED
LBL	LIGHT BLUE
BL	BLUE
LGR	LIGHT GREEN
GR	GREEN
VI	VILLET
OR	ORANGE
YE	YELLOW
WH	WHITE
LRD	LIGHT RED
PK	PINK
GY	GRAY
/	AND
I	OR

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER PROCESSES WRITTEN AUTHORIZATION
FROM VERMEIREN GROUP

EC A PLEGO VERSO MAIN PWR CABLE

Id: 1/1

Size: A4	Drafted	Verified	Comformed
Date:	2025-02-13	2025-02-13	2025-02-13
Initials:	FM	FM	TK

2998576/A

Plan No.

2998576/A

