


Use WEICONLOCK AN 302-43

Use WEICONLOCK AN 302-43

Use WEICONLOCK AN 302-43


Detail B
Scale: 3:5

Drawing shows only right side

BACKREST ADJUSTABLE ANGLE TRIGO				Plan No.
				EP00356/ J
Scale: 3:10	Drafted	Verified	Confirmed	
Date	29/10/2018	07/04/2023	07/04/2023	
Initials	PM	PM	TK	

EP00356 / J

Material				
No.	Plan	Article	Quantity	
* 1	YP06081(35)	2793265	1	
	YP06082(37)	2793365	1	
	YP06083(39)	2793465	1	
	YP06084(41)	2793565	1	
	YP06085(43)	2793665	1	
	YP06086(45)	2793765	1	
	YP06087(47)	2793865	1	
	YP06088(49)	2793965	1	
	YP06089(51)	2794065	1	
2	XP01632	GP09399	1	
3	XP03203	GP25700	1	
4	XP01862	GP13200	1	
5	M4x12	1907582	2	
6	---	---	---	
7	XP01861	GP12900	1	
8	M5x12 DIN 7991	2994076	2	
9	CAP Ø22	1990648	1	
** 10	Fold. push handle	EP00246	1	
** 11	M6x10	2994127	4	
	12	YP08333(35)	2613865	1
		YP08334(37)	2613965	1
		YP08335(39)	2614065	1
		YP08336(41)	2614165	1
		YP08337(43)	2614265	1
		YP08338(45)	2614365	1
		YP08339(47)	2614465	1
		YP08340(49)	2614565	1
		YP08341(51)	2614665	1
13	Hand grip GR-KS	2993809	1	
	GR-ERGO-S	2993810	1	
14	M6X6 DIN916/45H	1998138	1	



Plan	Lenght "X"
YP08333(35)	380mm
YP08334(37)	400mm
YP08335(39)	420mm
YP08336(41)	440mm
YP08337(43)	460mm
YP08338(45)	480mm
YP08339(47)	500mm
YP08340(49)	520mm
YP08341(51)	540mm

Plan	Lenght "X"
YP06081(35)	320mm
YP06082(37)	340mm
YP06083(39)	360mm
YP06084(41)	380mm
YP06085(43)	400mm
YP06086(45)	420mm
YP06087(47)	440mm
YP06088(49)	460mm
YP06089(51)	480mm

* When using a folded handle (10,11), select 2 sizes shorter tube.

** Option.

EP00356 / J

BACKREST ADJUSTABLE ANGLE TRIGO

Page 2/2

Plan No.

EP00356/ J

Scale: 3:10

Drafted

Verified

Confirmed

Date

29/10/2018

07/04/2023

07/04/2023

Initials

PM

PM

TK

